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► Se reporter au livret joint n°520937

[illegible]

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## CONTENTS

<b>1</b>	<b>INTRODUCTION .....</b>	<b>24</b>
1.1	Specification plates .....	24
1.2	Purpose of the device .....	25
1.3	Lifetime of the device .....	25
1.4	Target patient population .....	25
1.5	User group of the device .....	25
1.6	Limits of use .....	25
1.7	Side-effects on the patient .....	25
1.8	Warnings in case of device malfunction or variation in performance.....	25
<b>2</b>	<b>DESCRIPTION of an air plant HOSPITAIR: SCREW compressors.....</b>	<b>26</b>
2.1	Schematic diagrams .....	26
2.2	Description of the Air Production Assembly Architecture .....	27
2.3	Description of the Air Production Assembly Operation .....	27
2.4	Air treatment system.....	28
2.5	Receiver .....	28
<b>3</b>	<b>CHOICE OF THE PREMISES .....</b>	<b>29</b>
3.1	Environnemental conditions.....	29
3.2	Recommendation during installation.....	29
3.3	Installation .....	29
3.4	Ventilation of the premises .....	30
3.5	Fire.....	31
3.6	Packaging.....	31
3.7	Unloading and Installation .....	32
<b>4</b>	<b>INSTALLATIONS .....</b>	<b>32</b>
4.1	Dimensions .....	32
4.2	Electric power supply.....	34
4.3	Connections to the pneumatic network .....	35
4.4	Drains .....	36
<b>5</b>	<b>CONNECTIONS.....</b>	<b>36</b>
5.1	Installation of HOSPITAIR units, type G or Pack or Skid .....	36
5.2	Competencies, training and qualifications.....	36
5.3	Beginning of the work.....	37
5.4	Optional receiver .....	37
5.5	Pipes, flexible hoses and fittings .....	37
5.6	Brazing of piping .....	37
5.7	Checks during and after the installation of the pipes .....	38
5.8	Identification and traceability.....	38
<b>6</b>	<b>STARTUP PROCESS .....</b>	<b>38</b>
6.1	First step .....	38
6.2	Switching on .....	38

<b>7</b>	<b>SETTINGS.....</b>	<b>39</b>
7.1	Setting pressure .....	39
7.2	FLEXO time duration .....	39
7.3	Standby time duration .....	39
7.4	Hygrometry .....	39
<b>8</b>	<b>ALARMS .....</b>	<b>40</b>
<b>9</b>	<b>MAINTENANCE AND SERVICING .....</b>	<b>40</b>
9.1	Pressure indicator .....	41
9.2	Hygrometry probe .....	41
9.3	PAR filter replacement of final filtrations or pressure release line .....	41
9.4	Mechanical back-up service .....	42
9.5	Discharge system .....	42
9.6	Network alarm .....	42
9.7	Safety relief valve.....	42
9.8	PT100 temperature sensor .....	42
<b>10</b>	<b>24V OPTION with back-up for PROCOM 2.....</b>	<b>42</b>
10.1	Description.....	42
10.2	Start-up.....	43
10.3	Maintenance.....	43
<b>11</b>	<b>TEST PRESSURE .....</b>	<b>43</b>
<b>12</b>	<b>PROTECTION AGAINST OVERPRESSURE .....</b>	<b>43</b>
<b>13</b>	<b>TRAINING .....</b>	<b>43</b>
<b>14</b>	<b>CLEANING .....</b>	<b>43</b>
<b>15</b>	<b>END OF LIFE OF THE MEDICAL DEVICE .....</b>	<b>43</b>
<b>16</b>	<b>ELECTROMAGNETIC COMPATIBILITY .....</b>	<b>44</b>

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## MEANING OF THE PICTOGRAMS USED IN THIS MANUAL

Symbols are used in these instructions to help you understand them and to indicate particular points which you need to act upon or consider. Those symbols are:

**Reading symbol:** to refer to the operator manual / instructions manual.

**Danger symbol:** This symbol represents the DIN 4844 symbol (hazard warning); it warns of a hazard which could result in death or injury, along with damage to the machine. This symbol must be observed by all personnel working on this machine. Regulations on safety in the workplace must be observed.

**Environmental protection symbol:** This symbol serves as a reminder to sort waste during maintenance operations, to store it in a safe place, and to dispose of it with due care for the environment.

**Dustbin on wheels barred of a cross symbol:** This symbol indicates that the equipment should not be thrown with domestic waste and is subject to selective collection.

**Non-ionizing radiation hazard symbol:** This pictogram informs of a risk or a danger related to non-ionizing radiations.

**Medical Device:** The product is a medical device.

### General remark:

If used in accordance with their intended purpose, units bearing the CE mark will satisfy the essential requirements of directive 93/42/EEC, relative to medical devices and, where necessary, to essential requirements of the directive 2014/68/EU, relative to pressure equipment.

Those units, **when installed in accordance with these instructions**, comply with the provisions of currently applicable European Directives and harmonised standards.

However, your attention is drawn to the fact that a unit with only one compressor or only one vacuum pump cannot be used without an emergency backup device: gas backup or other fixed or mobile unit.

For devices intended to be fixed and connected to a medical distribution system, the CE marking provided by MIL'S as well as the CE conformity declaration of the product covered by the start-up and maintenance instructions becomes invalid for any relocation to another site/health establishment.

The device may only be relocated to another site if it has been reconditioned under the manufacturer's responsibility to comply with general safety and performance requirements, please contact us for more information.

To ensure personal safety and to avoid any damage to equipment, it is essential to follow the instructions contained in this document and any other documentation provided with the device, and particularly the instructions on "Safe Practice".

## 1 INTRODUCTION

We strongly advise you to read these installation instructions, before first switching on your air unit. The damage that could be caused by setting the machine up or installing it incorrectly would not be covered by our guarantee.

### 1.1 Specification plates

In order to simplify any future requests for spare parts or further information about your compressed air plant, we advise you to copy onto this page the specifications written on the plates fixed to the compressed air plant, as soon as you receive your machine.

Consult the copy of the specification plate fixed to the compressed air plant in fourth of cover.

Compressed air plant HOSPITAIR	
Type HOSPITAIR	
Model	
Regulation type	
Serial number	
Date of fabrication	
Power of the plant	
Power supply	
Maxi. Pressure	
Fluid group II	

Air treatment unit	
Type	SEC 7HC
Model	
Serial number	
Date of fabrication	
Power supply	
Maxi. Pressure	16 bar
Fluid group II	

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## 1.2 Purpose of the device

### 1.2.1 Medical indications

The device is intended to supply medicinal air as defined in European Monograph No. 1238. The air units are intended to produce continuously, within the pressure range required by the regulations, the air flow rate required for the care activities of all the patients in the healthcare establishment.

### 1.2.2 Medical contraindications

Not applicable.

## 1.3 Lifetime of the device

The device covered by the user manual has a lifetime of 10 years.

The performance of the device is maintained during its lifetime by:

- periodic verification and maintenance operations, in accordance with the recommendations of the circular of October 10, 1985,
- continuous analysis and measurement of performance with an alarm triggered in the event of a deterioration in performance,
- the replacement of defective components with original components that comply with the manufacturers' requirements,
- uncontrolled non-modification or non-extension of the distribution system,
- feedback on this type of equipment.

## 1.4 Target patient population

Patients do not physically handle and are not in direct contact with medical air units. These units are installed, maintained, and supervised by professionals. The medical AIR range is intended for all healthcare establishments requiring compressed air in their departments (breathing assistance, air instruments, drying of instruments after sterilisation, . . .).

## 1.5 User group of the device

The users of the "AIR" range depend on the life cycle of the product:

- **During the transport - storage phase**, the intended user profile is the transport operator. The device is packaged and therefore does not operate.
- **During the installation - connection, commissioning, shutdown, maintenance - servicing and dismantling phases**, the profile of the intended user is a qualified or technical department installer of the healthcare establishment. These are professionals who are familiar with the production and distribution systems for medical fluids.
- **During operation and use** the intended user is the technical department of the healthcare establishment. Similarly, these are professionals who are familiar with the production and distribution systems for medical fluids.

In order to ensure perfect aptitude in the use of the compressed air unit for medical use, the MIL'S SAS Company offers its own specific training for this equipment. It is intended for users and technical personnel working on or near the compressed air plant for medical use.

## 1.6 Limits of use

- See §.1.2.2.

## 1.7 Side-effects on the patient

The device has no clinical performance and therefore has no side-effects on the patient.

## 1.8 Warnings in case of device malfunction or variation in performance

The device has no clinical performance.

Medical personnel are not part of the user group claimed by the manufacturer. There is therefore no warning for medical personnel.

Cautions and warnings are indicated by the pictogram  in this manual.

In case of malfunction or variation in performance of the device, contact the manufacturer immediately.


→ Contact details on the cover page of this manual.

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type


## 2 DESCRIPTION of an air plant HOSPITAIR: SCREW compressors

### 2.1 Schematic diagrams

#### 2.1.1 11 bar version

#### 2.1.2 15 bar version

Note: the cyclone/filter and its electronic drain are included in the screw compressor starting from the MVD31 model.

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## 2.2 Description of the Air Production Assembly Architecture

The air is sucked in by the compressor through the air filter then is compressed and injected through the cyclone/filter fitted with an electronic air vent valve with an alarm contact. This one is used to trap the effluents before it enters the air treatment unit.

The pressure sensor of the compressor permanently controls the pressure before the air treatment unit.

Then the air enters the air treatment unit (see the description in the manual concerning the air treatment unit).

At the air treatment unit output, the air enters a manifold which supports the following items:

- ▶ PS : standby mechanical pressure switch of the regulation sensor,
- ▶ C1A : regulation sensor read by the PROCOM-2 or CYCLIC 2020,
- ▶ S1A : hygrometry probe.

For HOSPITAIR 15 bar version, the air is then stored in the receiver and then goes through the simple pressure release line [one pressure regulator + one particle filter PAR] or duplex parallel pressure release line [two pressure regulators + two particle filters PAR].

For HOSPITAIR 11 bar version, the air is then stored in the receiver and then goes through the simple final filtration [one particle filter PAR] or duplex parallel final filtration [two particle filters PAR]..

For duplex version, the air enters in two sets at the same time and each set can be "by-passed" for the maintenance operations with isolating valves.

- ▶ an analysis connector, to check the quality of the generated air.
- ▶ a network alarm sensor connector (C2A sensor).

## 2.3 Description of the Air Production Assembly Operation

Regulation is available either with PROCOM-2 or with CYCLIC 2020.

When the sensor C1A reaches 9 bar (11 bar version) / 10 bar (15 bar version), compressor A is forced into operation for a pre-set period of time. This forced operation time is called the FLEXO time.

The discharge system of the compressor is disabled, and the cycle of air dryer A is started.

### 2.3.1 During the FLEXO time:

The compressor is still operating, whatever the pressure at C1A. The pressure can exceed 9.5 bar (11 bar version) / 12.5 bar (15 bar version). If the pressure sensor of the compressor reaches 11 bar (11 bar version) / 14.5 bar (15 bar version), the compressor discharge system is enabled in order not to go over 11 bar (11 bar version) / 15 bar (15 bar version), then at 9.5 bar (11 bar version) / 13 bar (15 bar version), the discharge system is disabled again and the pressure starts rising again, and so on and so forth...

### 2.3.2 When the FLEXO time has elapsed:


Two situations are possible, as follows:

- a) The pressure at C1A level is up to 9.5 bar (11 bar version) / 12.5 bar (15 bar version):

The compressor is stopped according to its own procedure and an order is given to the air dryer to complete its cycle.

- b) The pressure at C1A level is between 9 and 9.5 (11 bar version) / between 10 and 12.5 bar (15 bar version):

When the pressure at C1A level reaches 9.5 bar (11 bar version) / 12.5 bar (15 bar version), the compressor is stopped according to its own procedure and an order is given to the air dryer to complete its cycle.

The next compressor to be started will be the B one, and so on and so forth.

### 2.3.3 Set-point values

#### Pd: pressure sensor of the compressor

- 11 bar version: is triggered at 11 bar and is reclosed at 9.5 bar.
- 15 bar version: is triggered at 14.5 bar and is reclosed at 13 bar.  
(connector on the screw compressor outlet / physical position : in the screw compressor, setting in the SIGMA CONTROL)

#### C1A: regulation sensor

- 11 bar version: 9 bar compressor start-up / 9.5 bar compressor stoppage after the FLEXO time.
- 15 bar version: 10 bar compressor start-up / 12.5 bar compressor stoppage after the FLEXO time.  
(connector on the receiver / physical position : regulation unit)

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## Ps: standby pressure switch

It is enabled to back-up the C1A sensor or the PROCOM-2 or CYCLIC 2020 board in case of a failure.

- 11 bar version: 7 bar starting of all compressors / 9 bar stoppage of all compressors.
- 15 bar version: 8 bar starting of all compressors / 10 bar stoppage of all compressors. (connector on the receiver / physical position: regulation unit)

## **2.3.4 Functionalities attached to the compressor**

A general alarm is transferred from the compressor to the PROCOM-2 or CYCLIC 2020.

This alarm is the synthesis of different status:

- ⇒ motor disjunction,
- ⇒ too high oil temperature,
- ⇒ compressor switched off,
- ⇒ wrong direction of rotation.

*Note: the compressor will start automatically after a breaking-off of circuit.*

## **2.3.5 Compressor emergency stop**

**Use the emergency button only in extreme necessity because it stops suddenly the compressor and if the compressor is on charge, oil can go back in suction.**

During test, it is better to use the I/O button of the SIGMA CONTROL to stop the compressor.

## **2.3.6 Information about shot risk**

The oil separator, mounting between gaskets, can be isolated from the ground of the unit. It can be load electrostatically and there is a risk of spontaneous ignition of the oil vapour. The oil separator gaskets are equipped with aluminium braid to maintain the ground. This situation is exceptional, but you must be informed to take safety precaution to operate and to use genuine material.

## **2.4 Air treatment system**

The air treatment system is composed of:

- A cyclone/filter fitted with an electronic air vent valve with an alarm contact in case of failure,
- A micronic filter, then a submicronic filter both equipped with condensate drain,
- An adsorption dryer with 3 levels of air treatment:
  - ⇒ Molecular sieve to remove CO<sub>2</sub> and humidity,
  - ⇒ Activated carbon to remove oil vapours and odours,
  - ⇒ Hopcalite to remove CO.
- an activated carbon filter for odour elimination.

Each air treatment system is designed to support the flow rate of one compressor.

## **2.5 Receiver**

**Service conditions:** The service conditions are specified on the receiver identification plate.

The operating pressure should never exceed the service pressure.

The following minimum and maximum temperatures shall be complied with:

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**These temperature values are given as a function of the capacity and only apply to the pressure receivers used with our HOSPITAIR products.**

**Purpose of use:** These are single pressure devices designed to contain air or nitrogen.

## **Installation and maintenance:**

The receiver must be installed in neutral, normal and current atmosphere and in places which are not likely to involve a rise in temperature of the device exceeding the temperature of service; it must be connected to ground. Access for check must be easy.

Incidents are always possible (exhaust of valve, leakage on drain or valve, etc...) so it's better to avoid the installation of a device in an attended site, or near fragile devices.

To check its interior state, the receiver must be visited regularly and according to the regulations in force in the country where it is installed.

A depth control must be practiced regularly with adequate means (ultrasound, magnetic...).

Fixing and fastening of receivers must be made while taking care to avoid all forced and tensions and to damp out all vibrations, even non-visible.

The base of the receiver will have to support the weight of the receiver, including in situation of tests.

In performing inspections, visits and periodical tests, we remind users to comply with the legislation in force where the tanks are in use in their country.

The equipment must not be modified without the manufacturer's prior approval.

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

**In complement of the indications hereafter, refer to the instruction manual provided with each receiver.**

The receivers and their accessories must be protected from freezing and the fire.

Use adequate handling and lifting system and take care not to damage the receiver or its components.

Before proceeding to the inspection and/or maintenance operations on the receiver and its components, make sure that this one is completely drained (internal pressure equal to the atmospheric pressure).

Not to carry out a welding or modifications (drilling, heating with a flame,...) on the parts under pressure of the receiver.

The connection of the receiver with the installation must be carried out with flexible connections in order to avoid all vibrations.

Take care that the receiver is always provided with its security accessories and its inspecting devices.

In the event of damage of the device, check by calling upon a specialist, that the receiver can be used in full safety.

## 3 CHOICE OF THE PREMISES

### 3.1 Environnemental conditions

**The ambient temperature should be between +5°C and +40°C.  
The room must be properly ventilated.  
The altitude must be between 0 and 500 meters.  
For operation out of theses ranges, please contact us as performance is impaired.**

Storage temperature: ..... +2°C à +60°C

Humidity: ..... 30 to 80% relative humidity

The medical air unit must be installed on a surface that can withstand a load drop of minimum 775 kg/m<sup>2</sup>.

The use of the compressed air units is exclusively in the premises of the healthcare establishment: technical room, room on the floor, outdoor cabin, . . . These premises are intended for the medical compressed air production. The medical air unit must be installed in a clean, non-humid, non-saline, dust-free environment.

### 3.2 Recommendation during installation

A back-up supply source composed of high-pressure cylinders (or cylinder racks) of medical air is recommended to ensure the continuity of the supply, even in the event of a major electrical failure.

The enclosures (interior or exterior) for the medical air feed systems must comply with the following requirements:

- When an enclosure is located near a source of heat such as an oven, incinerator or a boiler room, the construction must ensure that the temperature of the premises does not exceed 40°C.
- If a system for heating the premises is installed, no point in contact with the air must exceed 225°C.
- An enclosure must not be located within 3 m of an electric transformer,
- Enclosures must not be next to fuel storage tanks,
- Enclosures must have a concrete floor
- The fences and the walls must be at least 1.75m high
- The premises or zones occupied by the generators must not be used for other purposes.

The medical air production plant must be protected against shocks from passing vehicles, carts, etc.

Standards **NFS99-175** (§ 5-7-3) and **ISO 7396-1** (§ 5-8) indicate that an air plant must not be installed in the same room as the concentrators of O<sub>2</sub> or cylinders of non-cryogenic gases and liquids.

When choosing the location for a medical air production plant, it is necessary to take into account the potential risks (of contamination and fire for example). Drainage systems must be installed in this location.

The air admission(s) of the compressors must be in a place with minimum contamination by exhaust from an internal combustion engines (for example, automobiles), exhaust from vacuum systems, exhaust from medical gas distribution networks, from anesthesia gas exhaust systems, exhaust from ventilation systems and other sources of contamination.

The admissions must have means to prevent the entry of insects, debris and water.

It is important to take into account the potential effects of the prevailing winds on the location of the admission(s).

### 3.3 Installation

**Before any handling of the machines, see section §3.7.**

The installer bears the full responsibility for the environment where the product will be installed. Means for fighting fire, to protect the medical air production plant against shocks, four weather, floods, earthquakes and corrosion must be installed by the installer.

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type


All the constructions, whether buildings, civil engineering structures or earthworks are at the expense of the operator and/or of the installer.

The medical air production plant will only function properly on a horizontal surface that can bear the total mass of the system and its accessories. Be sure to install the medical air production plant and its accessories such that the nameplates are visible and accessible.

## 3.4 Ventilation of the premises

The medical air production plant must be installed in a sufficiently aerated room, free of dust and protected from frost. It can control a mechanical ventilation system slaved to the temperature.

### 3.4.1 Extraction fan

**A:** For the ventilation of the premises, we recommend renewal of 6 times the volume of the room per hour.

This flow  $Q_a$  should preferably be mechanically extracted.

$$Q_a = S \times h \times 6$$

With:

$Q_a$  = Air flow to be extracted (m<sup>3</sup>/h)

$S$  = Floor area of the premises (m<sup>2</sup>)

$h$  = Height of the premises (m)

**B :** If the cooling air discharge from the compressors is not channeled, this air flow to be extracted must be added according to the formula:

$$Q_b = \frac{P \times 1000 \times 3600}{\Delta T^\circ \times Cp \times \rho} = 3024 \times \frac{P}{\Delta T^\circ}$$

With:

$Q_b$  = Flow of air to be extracted (m<sup>3</sup>/h)

$P$  = Thermal power to be discharged (kW) – see table below.

$\Delta T^\circ$  = Temperature difference between the air in the premises and the outside air (10K maximum recommended)

$Cp$  = Specific heat of the air = 1005 J.kg<sup>-1</sup>.K<sup>-1</sup> at 25°C

$\rho$  = Density of the air = 1.1845 kg.m<sup>-3</sup> at 25°C

Type of compressor	Thermal power to be discharged by compressor (kW)	Heat generation (kCal/h)
MXV 2	2.2	1 892
MXV 3	3	2 580
MXV 4	4	3 440
MXV 5	5.5	4 730
MVA 9	7.5	6 450
MVB 12	11	8 340
MVB 16	15	10 576
MVC 19	18.5	18 920
MVC 23	22	18 920
MVD 31	30	25 800
MVE 38	37	31 820
MVE 46	45	38 700

**Example:** For two MVB12 compressors, with no sheathing of the cooling air discharge, in a 30 m<sup>2</sup> area with h=2.7 meters.

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

$Q_a$  : Renewal of the air of the premises 6 times per hour gives:

$$Q_a = 30 \times 2.7 \times 6 = 486 \text{ m}^3/\text{h}$$

$Q_b$ : Then, thermal power to be discharged:

$$P = 2 \times 11 \text{ kW} = 22 \text{ kW}$$

With allowable  $\Delta T = 10^\circ\text{C}$

$$Q_b = 3024 \times \frac{22}{10} = 6652 \text{ m}^3/\text{h}$$

The sum of the 2 flow is:  $6652 + 486 = 7138 \text{ m}^3/\text{h}$ .

The ventilation must be at least **7138 m<sup>3</sup>/h**.

## 3.4.2 Air intake

For the dimensioning of the air intake areas, the passage speed must not exceed 3 m/s and three flows must be included ( $Q_a + [Q_b \text{ or } Q_c] + Q_d$ ):

$Q_a$  : The flow of renewal of the surrounding air  $Q_a$  calculated as in the example of point 3.5.1.

$Q_b$  : If the cooling air discharge from the compressors **is not channeled**, you must add the air flow  $Q_b$  calculated as in the example of point 3.5.1.:

or

$Q_c$  : If the cooling air discharge from the compressors **is channeled**, the flow aspirated by the cooling of each compressor according to the table below must be added

$Q_d$  : The consumption of air aspirated by each compressor according to the table below:

Type of compressor	Thermal power to be discharged by compressor (kW)	Cooling air flow from the compressor (m <sup>3</sup> /h)	Air flow aspirated by the compressor at 11 bar (m <sup>3</sup> /h)	Air flow aspirated by the compressor at 15 bar (m <sup>3</sup> /h)
MXV 2	2.2	1000	16	
MXV 3	3	1000	22	15
MXV 4	4	1000	29	22
MXV 5	5.5	1300	40	32
MVA 9	7.5	2100	60	45
MVB 12	11	2500	100	78
MVB 16	15	3000	126	102
MVC 19	18.5	4000	179	150
MVC 23	22	5000	210	174
MVD 31	30	5400	269	223
MVE 38	37	8000	336	266
MVE 46	45	8000	411	328

The air intake area to be respected is given by the following formula:

$$S = \frac{Q}{V \times 3600}$$

With:

$S$  = air intake area (m<sup>2</sup>)

$Q$  = sum of the ventilation air flows (m<sup>3</sup>/h)

$V$  = passage speed (m/s)

**Example:** For 2 MVB12 – 15 bar compressors, with no sheathing of the cooling air discharge, in a 30 m<sup>2</sup> area with h=2.7 meters.

$Q_a$  : Renewal of the air of the premises 6 times per hour gives:

$$Q_a = 30 \times 2.7 \times 6 = 486 \text{ m}^3/\text{h}$$

$Q_b$ : Then, thermal power to be discharged:

$$P = 2 \times 11 \text{ kW} = 22 \text{ kW}$$

With allowable  $\Delta T = 10^\circ\text{C}$

$$Q_b = 3024 \times \frac{22}{10} = 6652 \text{ m}^3/\text{h}$$

$Q_c$  : Not applicable

$Q_d$  : The consumption of air aspirated by each compressor

$$Q_d = 78 \times 2 = 156 \text{ m}^3/\text{h}$$

Whence:

$$S = \frac{Q}{V \times 3600} = \frac{486 + 6652 + 156}{3 \times 3600} = 0,67 \text{ m}^2$$

In this example, the real area of passage of the intake of air from the premises must therefore be at least 0.67m<sup>2</sup>, not including the taking into account of filters or anti-insect screens.

## 3.5 Fire

Suitable firefighting equipment must be installed and maintained in good condition in the premises of the generators and of the linked equipment in order to limit in particular the risk linked to the tanks of oxygen under pressure.

## 3.6 Packaging

The machines must be correctly protected in their original packaging, before being transported towards their final destination.

**Be careful using the cutter, because there is a risk of damaging the screen of the PLC, the paint, the tubes or the cables. In the event of deterioration, please contact the manufacturer.**

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## 3.7 Unloading and Installation

### 3.7.1 Compressors / Dryers

- ▶ Refer to the notices for the compressors and the dryers.

**Be careful when handling on an inclined plane.  
Handling with the help of a suitable pallet truck can require the strength of several people, even on a horizontal surface.**

### 3.7.2 Receivers

Each air receiver must be installed in a neutral atmosphere and in premises which would not cause an overall temperature increase of the unit exceeding the service temperature. It must be grounded.

Access for verification must be easy.

As incidents can always happen (leak on pipe or valve, etc.) it is best to avoid installing a unit in a very busy place or near fragile equipment.

The fastening and flanging of the receivers must be done making sure to avoid all stresses and tensions and to damp all vibrations, even those that are not apparent.

The base of the receiver must bear the weight of the tank, even in re-testing situations.

The receiver must not be modified or repaired without the manufacturer's prior authorization.

## 4 INSTALLATIONS

The installer and his employees must be able to prove that they have good knowledge of medical gases, their use and the associated risks. This level of knowledge can be demonstrated by initial training, additional training programs and experience.

The personnel working on this project must be able to demonstrate, always, that they have the required regulatory qualifications and certifications (electrical certification in particular).

### 4.1 Dimensions


#### 4.1.1 Compressors

- ▶ To refer to its technical manual attached.

#### 4.1.2 Air treatment unit SEC 7HC – AD5XXX

- ▶ To refer to its technical manual attached.


#### 4.1.3 Vertical air receivers


	11 bar			15 bar		
Capacity (litres)	500	900	1500	500	1000	1500
Service pressure	11 bar	11 bar	12 bar	15 bar	15 bar	15 bar
Ø	600	795	1000	650	800	950
A	2081	2153	2305	1850	2350	2550
B	1655	1798	1780	1463	1756	1660
C	1255	1370	1468	1243	1556	1460
E	785	873	680	443	506	560
F	173	153	170	123	96	80
Inlet valve	1"	1"	1"1/4	1"1/4	1"1/4	1"1/4
Outlet valve	1"	1"	1"1/4	1"1/4	1"1/4	1"1/4

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## 4.1.4 Simple or duplex parallel pressure release lines for 15 bar version

Dimensions in mm

Note: the standard assembly in wall fixation version allows a passage of the medical air from left to right. The direction can be easily reversed by turning the brackets.

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

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## 4.2 Electric power supply

See the electric diagram attached to this manual (terminals and electrical connections).

All operations on the air plant shall be performed with the air plant switched off.  
Connection must be carried out by the operator.

Connect at least one air production assembly to a backed-up electric power source, as well as the regulation unit.  
If two backed-up electric power sources are available, it is recommended to connect the various air production assemblies to different power sources.

Class I device: to avoid any risk of electric shock, this device must only be connected to a power supply network equipped with a protective earth.

Check that the cables and wires inside the electric cabinet are clamped tight to their terminals before the initial start-up.

- Check that the voltage read on the motor plates or the station plate matches the voltage of your power system.

Protect all power leads located upstream by an over-current protection device in compliance with the legislation and regulations in force in your country.

The wire sections are given for information.

Even so, they shall be defined by taking contacts with the wire provider.

They depend:

- ⇒ of the wire type
- ⇒ of the length
- ⇒ of the room temperature
- ⇒ of the wire path

### 4.2.1 Compressor power supply

- Each compressor should be fed independently, in terms of power, by 3 wires + ground (from the general electric supply system).

### Protections in Ampers

Screw compressor 3-ph 400V 50Hz	MX 2	MX 3	MX 4	MX 5	MVA 9	MVB 12	MVB 16	MVC 19	MVC 23	MVD 31	MVE 38	MVE 46
Power (kW) for each compressor	2.2	3	4	5.5	7.5	11	15	18.5	22	30	37	45
Strength of current (A)	4.9	6.4	8.1	10.8	15	24	30	40	46	62	74	91
Wire cross section (mm <sup>2</sup> ) for each compressor	2.5	2.5	2.5	2.5	4	6	10	10	16	25	35	35
Protection to install * (A)	10	16	16	16	25	32	40	63	63	80	100	100

\* A D-curve circuit breaker protection must be installed upstream the compressor.

### 4.2.2 PROCOM-2 or CYCLIC 2020 power supply

Feed the PROCOM-2 or C2020 electrical box with 1-ph 230V – 50/60Hz.

If only one compressor is present, a 2 x 0.75 mm<sup>2</sup> + earth cable is sufficient.

If at least two compressors are installed, provide two 2 x 0.75 mm<sup>2</sup> + earth cables. In this case the power supply X must come from the same source as that of compressor A and the power supply Y from the same source as that of compressor B.

*Note: Connect the PROCOM-2 or CYCLIC 2020 alarm synthesis (AL1 / AL2) to a visual and/or audible alarm on the control unit.*

### 4.2.3 Connection between compressor and PROCOM-2 or CYCLIC 2020 electrical box (only HOSPITAIR PACK piloted and HOSPITAIR G)

Each compressor is delivered with a 10-meter long cable with a screw connector to connect to the PROCOM-2 or CYCLIC 2020 electrical box.

The following information is provided:

- ⇒ The PROCOM-2 sends information to the compressor to start and to stop,
- ⇒ The compressor sends compressor defect or Beko drain defect information to the PROCOM 2.

Screw the connector on the plug underneath the regulation unit PROCOM-2.

*Note: if the 10-meter long cable is not enough, an extension can be added. Contact us.*

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type


## 4.2.4 Connection between air treatment system and compressor (only HOSPITAIR G)

Each air treatment system is delivered with two 10-meter long cables with a screw connector to connect to the compressor.


- First connector M12 to feed in 230V the dryer and the Beko drain from the compressor.
- Second connector M8: the compressor sends information to the dryer to start and to stop and the Beko drain sends information of defect to the compressor.

## 4.2.5 Summary of electric connections

### HOSPITAIR PACK et SKID

### HOSPITAIR G

## 4.3 Connections to the pneumatic network

The connections located after the air treatment units are essential. These links must be visited and inspected regularly.

After the pneumatic connection of all the elements between them, check the sealing of each flexible connections.

### 4.3.1 Connection compressor / air treatment system (only HOSPITAIR G)

- Connect the compressor outlet to the air treatment unit, using the pressure hose.

### 4.3.2 Connection air treatment system / air receiver / final filtration or pressure release line

- Connect each air treatment unit outlet to the final filtration or pressure release line on the network.
- Connect air receiver(s) in dead end between the air treatment system outlet and the final filtration or pressure release line inlet.

*Note: it is possible to connect the air receiver in crossing position (in line) between the air treatment unit and the final filtration or pressure release line. The crossing position makes the average of CO2 rate possible, but a by-pass system is needed.*

### 4.3.3 Control unit

- Connect on the left side of the PROCOM 2 or CYCLIC 2020 electrical box, where there are 3 sliding valves, the PTFE white or copper or stainless-steel pipe dia.6 from the tapping connection between the air treatment units and the final filtration or pressure release line to deliver air:
  - ⇒ To the regulation sensor C1A,
  - ⇒ To the emergency mechanical pressure switch,
  - ⇒ To the hygrometry sensor S1A.

### 4.3.4 Network alarm

The network sensor C2A is fitted in the PROCOM 2 or CYCLIC 2020 electrical box.

- Connect on the right side of the PROCOM 2 or CYCLIC 2020 electrical box, where there is a sliding valve, the RILSAN black or copper or stainless-steel pipe dia.4 from the final filtration or pressure release line to deliver air to the network sensor C2A.

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## Alarm report:

The connection between Alarm report on automate units must be done by the person in charge of the Installation.

He must control the compliance wiring and proper operation during the final acceptance test.

## 4.4 Drains

The effluents coming from the PF and SUB filters' drain valves and the electronic vent valve located in the compressor must be collected.

In order to respect the environment, these condensates must be treated before discharge to the sewer. (oil and water phase separation by using an OWAMAT separator)

In accordance with:

Law 76-663 of 19 July 1976 on installations classified for environmental protection.

Law 92-3 of 3 January 1992 on water.

Decree 93-742 of 29/03/93 on the authorisation and declaration procedures provided for in Article 10 of Law 92-3 of 3 January 1992.

Decree 93-743 of 29/03/93 on the nomenclature of transactions subject to authorisation or declaration pursuant to Article 10 of Law 92-3 of 3 January 1992 on water.

Order of 1 March 1993 on withdrawals and consumption of water and discharges of all kinds from classified installations for the protection of the environment subject to authorisation.


Note: for installing the OWAMAT separator (option), refer to the OWAMAT installation and operation manual attached hereto, and strictly comply with all recommendations.

## 5 CONNECTIONS

### 5.1 Installation of HOSPITAIR units, type G or Pack or Skid

The HOSPITAIR compressed air unit packages in G or Pack or Skid configuration are delivered as separate elements which must be connected on site by the installer or the customer.

The diagram below shows the connection operations to be carried out:

The connecting hoses between the compressors and the dryers are provided by the MIL'S SAS Company.

The connections (red link in the diagram) between the dryers, the final filtration line and the optional tank are at the expense of the installer or the client, as is the electrical connection of the various components.

The following sections detail the requirements taken from the risk analysis which must be observed by the installer or the client during installation and connection:

### 5.2 Competencies, training and qualifications

The installer and his employees must be able to prove that they have good knowledge of medical gases, their use and the associated risks. This level of knowledge can be demonstrated by initial education, additional training programs and experience.

The personnel working on this project must be able to demonstrate, always, that they have the required regulatory qualifications and certifications (electrical certification in particular).

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

It is imperative that the installer:

- qualify his personnel who do brazing in compliance with the requirements of the **ISO 13585** standard,
- qualify his brazing procedures in compliance with standard **EN 13134**,
- keep records concerning the qualifications of the personnel and procedures for a period of at least 15 years. These records must be presented at the request of the MIL'S Company, the Competent Authorities, and the notified body of the manufacturer.

## 5.3 Beginning of the work

Work may not begin without the issuing of a Work Authorization from the operator.

If a network shutdown is necessary, the sub-contractor MUST follow the shutdown procedure in effect in the establishment. The shutdown cannot be carried out unless the operator has issued a shutdown authorization.

The sub-contractor must comply with the requirements of the health care establishment in terms of hygiene and cleanliness, particularly during work in sterile areas and/or areas occupied by patients. Special attention will be paid to the risk of bio-contamination of equipment and components.

Before the installation of each component, the sub-contractor must verify its state of integrity, its level of cleanliness and, if applicable, the best by or expiry date. All components which have damaged packaging must be discarded.

The sub-contractor must be particularly vigilant, during work near sources of production, regarding the possible automatic restarting of rotating machines (remote-controlled starting).

For work on electrical control panels, lock-out systems must be used.

All of the components used must be compatible with oxygen and delivered clean, in compliance with the requirements of standard **NF EN ISO 15001**.

## 5.4 Optional receiver

If the operator takes charge of supplying the storage receiver, it must be adapted to the nominal working pressure of the facility and must bear the CE marking in compliance with directive **2014/68/UE** (equipment under pressure) or **2014/29/UE** (simple pressure vessels).

The receiver must be equipped with:

- one or several sectional valves,
- a CE marked safety valve in compliance with directive **2014/68/UE** (equipment under pressure),
- a manometer,
- a drain valve.

## 5.5 Pipes, flexible hoses and fittings

The rigid pipes must be made of copper in compliance with standard **NF EN 13348**.

The pipe diameters must be identical to the input / outlet diameters of the elements to be connected. The piping, of a diameter less than 32 mm (diameter 1" ¼ maximum) is, in the meaning of directive 2014/68/UE, classified in the meaning of article 4.3.

The fittings used to connect the copper pipes must be compliant with **EN 1254-1** and **EN 1254-4**.

If low pressure connecting hoses are used, they must be in compliance with standard **NF EN ISO 5359**.

All of the components used must be delivered clean, protected from all contamination before and during the installation.

Furthermore, the components and materials must be able to withstand 1.5 times the operating pressure for the compressed gases.

Proof of the maintaining of the pressure of the piping and components used must be presented by the installer upon request from the MIL'S Company, the Competent Authorities or the notified body of the manufacturer.

The storage and handling conditions must ensure the preservation of the product and of its level of cleanliness, until the final assembly.

The cutting of piping intended for the distribution of compressed medical fluids must be **EXCLUSIVELY** done with a pipe cutter in order to limit the risk of introduction of shavings or particles.

The installation of the piping must be done according to Professional Standards and must meet the applicable requirements of section **§11.1** of standard **NF EN ISO 7396-1**.

The pipe supports must be compliant with the requirements of section **§11.2** of standard **NF EN ISO 7396-1**.

All of the components of the network must be connected to a ground connector in compliance with the regulations in effect.

The pipes must be protected against all physical damage.

## 5.6 Brazing of piping

With the exception of the mechanical fittings used for certain components, all of the metallic fittings must be brazed or welded. If a filler metal is used, its temperature in the liquid state must not be less than 600 °C. The filler metals must be nominally free of cadmium. If silver alloys are used, they must be compliant with **ISO 17672**.

During the brazing or welding of the pipe fittings, the interior of the pipe must be constantly purged with a protective gas in order to limit the risk of introduction of shavings or particles.

The sub-contractor cannot carry out brazing operations unless he has a fire permit issued by the health care establishment.

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## 5.7 Checks during and after the installation of the pipes

Self-checks must be carried out during the installation, particularly a visual verification of the brazing in order to detect:

- the absence of flux during the brazing operation,
- any deterioration of the base metal due to overheating,
- the absence of filler metal at the joint,
- the absence of run-out or surplus filler metal,

The **network piping and components providing the link between the various components** of the production unit must be subjected to a tightness test at 1.43 times the maximum allowable pressure. The table below indicates the pressure to be applied for this test:

Pressure (bar)	Test pressure (bar)
7	10.5
14	21

The results of the verifications must be recorded and saved.

The entire facility must undergo a final verification before commissioning. Standard **NF EN ISO 7396-1** applicable to medical gas distribution systems details the tests to be carried out.

The results of the verifications must be recorded and kept.

## 5.8 Identification and traceability

There must be traceability for the following elements.

- Sources of production,
- Storage receiver,
- Safety valves installed on the network,
- Valves (if they bear a serial or batch number).

The sub-contractor **MUST** indicate the serial or batch numbers of the above components on the self-verification report.

The sub-contractor must label the distribution system and the valves at risk in compliance with the standards and regulations (**ruling of December 10, 2004** in particular) in effect.

## 6 STARTUP PROCESS


**Screw compressor which will be stored for more than 6 months must be drained; air and oil filters and oil separator must be replaced.**

In the contrary case, the filters tend to dry and create an abnormal pressure loss causing a pressure higher than the valve opening pressure of the separator tank of the compressor and a very strong oil mist can escape. Oil can degrade and more not provide its function.

## 6.1 First step

- Check the oil level of compressors.
- Check that the button AUTO/CONTINU in the front side of the compressor is in the AUTO position.

*Note: the button on the CONTINU position allows each line to operate in forced running.*

**If the button of a compressor is set to CONTINUOUS, non-conforming air can be sent to the network. The H<sub>2</sub>O/CO/CO<sub>2</sub>/O<sub>2</sub>/NO/NO<sub>2</sub>/SO<sub>2</sub> gas analysis alarms will not be able to stop the active compressors, despite the effect setting to STOP in the PLC.**

- Before starting the equipment, first read all the running operation and maintenance manuals (screw compressor, ADSORMIL'S air dryer, PROCOM-2, CYCLIC 2020).
- Check, with the electrical drawing, the switches of PROCOM-2 or CYCLIC 2020 (right configuration).
- Set the On/Off switch to O (Off) on all compressors.
- Press simultaneously the  and  keys of the CYCLIC 2020 to mask the hygrometry during 18 hours. The °C LED flashes during this time.
- Hygrometry is masked during 24 hours on the PROCOM-2.

## 6.2 Switching on

*Warning:*

*Upon first start-up, or if a major incident has occurred on the network and the network had to be brought back to atmospheric pressure, during the first moments of operation, the quality of the air will not be in compliance with the european monograph (past: NFS 90 155 standard), as the ambient air which was present in the whole system has to be eliminated. Therefore the pressure reducing valve outlet must be shut off, pressure must be increased, the compressors must be stopped, and the system must be vented until a residual pressure of 2 bar is obtained in the receiver (this procedure is described below).*

*Note: during the first commissioning procedure, or after an intervention on the system, these operations will not solve the hygrometry problem.*

Each screw compressor is equipped with a phase controller indicates phases are in correct position. If not, the compressor can't start.

- Shut off the network outlet valve on the final filtration or pressure release line outlet.

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

- Start the compressor A pushing the "I" key on the SIGMA CONTROL: if there is an electrical problem and fault, the compressor stops immediately. In this case, insulate the compressor electrically and cross 2 wires at compressor motor level
- Stop the compressor A and start compressors B and C to check the rotation direction of each compressor.
- Restart all compressors to allow the pressure to go up to maximum level.
- Check the links for leaks.
- For the 15 bar version, check the pressure obtained downstream of the pressure release line outlet, create a consumption through the analysis valve and use the reducing valve to adjust the pressure as appropriate (pressure setting: 9.5 bar).
- Cut off all compressors and vent the receiver to 2 bar.
- Restart all the compressors and open the network outlet valve gradually.
- To check the regulation for correct operation:
  - ⇒ isolate the system from the network, using the final filtration or pressure reducing valve,
  - ⇒ Create a leak that you will control according to the intended purpose, on the receiver vent valve,
  - ⇒ After the check, shut off the vent valve and open the network valve.
- For further information concerning the PROCOM-2 or CYCLIC 2020 regulation mode, refer to their manuals.

## 7 SETTINGS

Any modification to the settings should only be carried out by qualified and authorised personnel.  
Improper actions may cause complete system malfunction.

### 7.1 Setting pressure

#### 7.1.1 For CYCLIC 2020 / PROCOM-2: C1A

Factory setting	11 bar Version	15 bar Version	Parameters CYCLIC 2020
Upper limit	9.5 bar	11.5 bar	Parameter 0
Lower limit 1 <sup>st</sup> compressor	9 bar	10 bar	Parameter 1
Lower limit 2 <sup>nd</sup> compressor	8.5 bar	9.5 bar	Parameter 2
Lower limit 3 <sup>rd</sup> compressor	8 bar	9 bar	Parameter 3
Regulation alarm	7.5 bar	8.5 bar	Parameter 4

#### For CYCLIC 2020 :

**Modification:** Press keys  and  simultaneously for 2 seconds. Use the key  or  to select parameter 0, 1, 2, 3 or 4 and adjust the value using the key  or .

*Note: repeat the operation for the other parameters.*

**For PROCOM-2 :** modify values in the parameter page.

#### 7.1.2 Standby mechanical pressure switch: Ps

Factory settings	11 bar Version	15 bar Version
Maxi	9 bar	10 bar
Mini	7 bar	7 bar

#### Modification:

Unscrew the green screw to move the green indicator down.  
To set the stop threshold: adjust the red indicator near 10 bar, screw (if to low) or unscrew (if to high) the red screw.  
To set the hysteresis of 2 bar, screw or unscrew the green screw.

### 7.2 FLEXO time duration

**Factory setting:** PROCOM-2: 10 min / CYCLIC 2020: 8 min.

**Modification:** use switches SW2-1 and SW2-2 (for CYCLIC 2020) or act in the parameters page (for PROCOM-2).

### 7.3 Standby time duration

**Factory setting:** 30 minutes (G2/G3) – unlimited (G1).

**Modification:** use switches SW2-4 and SW2-5 (for CYCLIC 2020) or act in the parameters page (for PROCOM-2).

### 7.4 Hygrometry

**Factory setting:** -45°C (parameter 8).

**Modification:**

**For CYCLIC 2020:** Press keys  and  simultaneously for 2 seconds. Use the key  or  to select parameter 8 and adjust the value using the key  or .

**For PROCOM-2:** modify values in the parameter page.

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## 8 ALARMS

### For PROCOM-2:

The alarm page displays the alarms that are still present on the system.

### For CYCLIC 2020:

When the LED on the  key is lit, an alarm has occurred.

	constant	fault on concerned compressor, may be in connection with: ⇒ Disjunction ⇒ Drain under the cyclone			
	flashing	hygrometry fault on concerned line			
	flashing	upper limit fault on concerned compressor			
	constant	temperature > 40°C in the room			
	Flashing slowly	temperature < 5°C in the room			
	constant	hygrometry fault			
			Version	11 bar	15 bar
	flashing slowly	regulation pressure fault	(parameter 4)	7.5 bar	8.5 bar
	flashing rapidly	back-up cylinder half empty or cylinders activation	(parameter 5) (parameter 9)	100 bar 7 bar	100 bar 7.5 bar
	constant	mini network pressure fault maxi network pressure fault	(parameter 6) (parameter 7)	4.5 bar 10 bar	6.3 bar 12.5 bar

Acknowledgement: press  and  simultaneously to acknowledge the alarm once read.

## 9 MAINTENANCE AND SERVICING

### IMPORTANT:

For this operation, use all necessary protective measures (gloves, mask, overalls etc.) then destroy contaminated consumables (antibacterial cartridges, condensates etc.) in according with that country's applicable laws (regulation incinerator).

Maintenance is to be entrusted to trained and approved personnel by the exploiting, with all necessary knowledge and acquainted with best practice.

Not to carry out maintenance operations when the device is in use.

If the spare parts used during maintenance operations are not corresponding to the manufacturers' requirements, the manufacturer declines all responsibility in the event of an incident.

The modification of the compressed air plant is forbidden

Modification of the compressed air plant is forbidden without the approval of the manufacturer.

If the compressed air plant is modified following manufacturer agreement, a control and a suitable test must be realized to make sure that the compressed air plant is always safe to use.

To avoid risk of damage or prevent equipment become dangerous is essential use the compatible components provide by the manufacturer.

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

Main components	6 months	12 months	24 months	3 years	5 years
COMPRESSOR		Refer to the maintenance manual			
AIR TREATMENT SYSTEM		Refer to the maintenance manual			
- Sequential / electronic drain test	•				
PRESSURE INDICATORS (sensor, pressure gauge)					
- Measurement validity check			•		
- Pressure gauge change					•
HYGROMETRY PROBE					
- Measurement check or re-calibration				•	
AIR RECEIVER					
- Call upon an inspection agency		As per the legislation in force			
FINAL FILTRATION					
- Particles filter (PAR) change			•		
PRESSURE RELEASE LINE					
- Particles filter (PAR) change			•		
- Pressure reducing valve check		•			
- Regulator maintenance				•	
ELECTRICAL BOX					
- Clean and dust using a soft brush		•			
- Check the connections in the box		•			
- Check the connections at motor terminals		•			
- Check the intensity of the motors		•			
- Operate the mechanical back-up system		•			
- Check the discharge system		•			
NETWORK ALARM					
- Check the general condition		•			
HOSE OR RILSAN UNDER PRESSURE					
- Overall condition check		•			
- Change					•
SAFETY VALVE					
- Change					•
DRAIN HOSE OR RILSAN					
- Overall condition check		•			
- Change			•		
OWAMAT (see the attached OWAMAT manual)		Replace the cartridge when water becomes unclear			

Note: The maintenance periodicities above are given for reference only and may be increased as required.

Note: Maintenance tables are available for the compressor, the air treatment unit and the pressure reducing valve outlet.

## 9.1 Pressure indicator

### Sensor


Connect a master pressure gauge on the test valve and open it.  
Check the measurement.

Note: the exhaust sliding valve allows:


- the sensor to be put in communication with the network,
- a pressure drop to be generated, in the form of a small leak, on the sensor.

During operation, the sliding valve must always be in full-passage position with the network.

normal position

test position (intermediate position)

**Pressure gauge:** a monitoring connector is fitted near every pressure gauge.

## 9.2 Hygrometry probe

The hygrometry probes slightly drift in time. Therefore they must be compared to a reference probe in our factory or sent back for re-calibration.

## 9.3 PAR filter replacement of final filtrations or pressure release line

### 9.3.1 Simple final filtration or pressure release line

Isolate filter with upstream and downstream valves.

Depressurise this part with vent valve.

Replace the cartridge. Re-install the assembly and set the valves in filtration position.

This operation will cause a rupture of the production. To make sure that another source will make it possible to ensure the continuity of production during this operation.

### 9.3.2 Duplex parallel final filtration or pressure release line

Isolate filters with upstream and downstream valves.

Depressurise this part with vent valve.

Replace the cartridges. Re-install the assembly and set the valves in filtration position.

Repeat this operation for the other part.

### 9.3.3 Checking the pressure regulator of the pressure release line

Setting pressure: 9.5 bar.

Create a leak downstream the pressure reducing valve and check the stability of it.

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## 9.4 Mechanical back-up service

Ask the cognizant authority whether this test is required or not, according to the activity and architecture of the facilities concerned. This action will trigger an alarm.

A mechanical standby pressure switch is mounted in parallel with sensor C1A on the PROCOM-2 or CYCLIC 2020.

- ⇒ set the master switch to "O".
- ⇒ remove the PROCOM-2 or CYCLIC 2020 power supply terminal box to cut the power supply.
- ⇒ set the master switch back to its initial position.
- ⇒ reduce the pressure. The compressors will start running at ~ 7 bar (11 bar version) / 8 bar (15 bar version) and will stop at around 9 bar (11 bar version) / 10 bar (15 bar version).

After performing this check, set the master switch to "O", put the PROCOM-2 or CYCLIC 2020 power supply back into place and reset the master switch to "I".

## 9.5 Discharge system

When the compressor rotates, make sure that the compressor discharge system operates when pressure reaches 11 bar (11 bar version) / 15 bar (15 bar version) or when an order is given to the compressor to stop, the compressor discharge system is activated.

Check the compressor discharge time before it stops completely (10 to 15 seconds).

## 9.6 Network alarm

Ask the cognizant authority whether this test is required or not, according to the activity and architecture of the facilities concerned. This action will trigger an alarm.

Each network alarm is fitted with an exhaust sliding valve.

To test the alarm, take the sliding part away from the nut part of the valve: a leak is created at sensor level; the alarm is triggered. Upon completion of the test, put the sliding valve back in normal position.

## 9.7 Safety relief valve

When a relief valve with discharge line operates and must be tested, no persons may stand in front of the valve discharge direction. Not allow discharge toxic, explosive or flammable products, into the atmosphere.

Before the test, a controlled discharge procedure in a confined volume is requested.

Never tamper with the valve, nor remove the lead/manufacture's seal, nor modify its adjustment pressure for any reason.

Not create cold or hot thermal shock on the relief valve.

Make sure that the lead/manufacture's seal is never damaged.

Maintain the valve clean and in functioning order. Check in particular that the exhaust remains free and that no foreign body blocks the evacuation piping.

To make sure that the safety valves continue to operate in good working condition, these must be periodically tested.

To do this, open the valve manually using the lever or the test knob; this test must be done while keeping the protected appliance at a pressure between 80 and 90% of the valve calibration pressure. The valve must open cleanly, and release an abundant amount of fluid, and must then close fully once the lever has been released or the knob tightened.

## 9.8 PT100 temperature sensor

It is recommended to change the PT100 temperature sensors every 5 years.

## 10 24V OPTION with back-up for PROCOM 2

### 10.1 Description

The 24V DC backup power supply module comprises an ELC backup module, two 12V DC batteries assembled in series to produce 24V DC and a 24V DC power supply adjusted to deliver a voltage of 27.3V DC.

Its purpose is to eliminate micro-power cuts on the 24V DC power supply to the device.

The purpose of the ELC backup module is:

- ⇒ to charge the batteries, controlling the charging current (0.3A)
- ⇒ to switch the 24V DC power supply over to battery power in the event of a 230V power supply failure.
- ⇒ to check that the batteries are connected the right way round

# HOSPITAIR PACK, SKID and G with SCREW (MV) – 11 bar and 15 bar PROCOM 2 or CYCLIC 2020 regulation - SEC 7HC air treatment type

## 10.2 Start-up

Switching on the power is sufficient. The ELC module will control the charging of the batteries.  
Charging time can be up to 8 hours.

Checking proper operation:

- ⇒ On receipt of the machine and after charging the batteries, disconnect the 230V power supply protection fuse of the ELC module and check that the PROCOM 2 is still live.
- ⇒ Every year, check that the backup power supply is in good working order using the procedure above.

## 10.3 Maintenance

The estimated length of service of the 2 batteries is at least 5 years. They must be replaced after 4 years in order to ensure that the backup power supply works properly.

Spare battery kit: **625144**

**Maintenance work should never be carried out when the power is on.  
Maintenance must be carried out by a qualified electrician.**

For this purpose:

- ⇒ Isolate the circuit electrically
- ⇒ Unplug the electric cables, locating the + and the -
- ⇒ Unscrew the battery support plate and replace it with a new assembly
- ⇒ Plug the electric cables back in, ensuring the right battery polarity
- ⇒ Turn the power back on and check for proper operation.

**Dispose of the batteries in accordance with the regulations in the country.**

## 11 TEST PRESSURE

The pressure equipment has been individually tested by its respective manufacturer.  
The equipment is connected to each other by means of a welding-free connection.

## 12 PROTECTION AGAINST OVERPRESSURE

The various modules constituting the compressed air unit, object of this technical manual, are protected by safety relief valves as indicated in the table hereafter:

Component	11 bar version Pressure setting (bar)	15 bar version Pressure setting (bar)	Location of the safety relief valve
Compressor	13 bar	16 bar	In the compressor
Dryer	---	---	Protected by the safety relief valve of the compressor
Receiver	11 bar	15 bar	On the receiver
Pressure release line	---	14 bar	After each regulator

## 13 TRAINING

In order to ensure perfect aptitude in the use of the compressed air unit, the MIL'S Company offers its own training for the use of this equipment. It is intended for users and technical personnel working on or near the compressed air unit.

## 14 CLEANING

**Do not use corrosive product.**

If the unit system requires cleaning, use a soap solution and a clean cloth.  
Then wipe well with a dry cloth.  
For disinfection, use swabs soaked in a disinfectant solution.  
Frequency: as required.

## 15 END OF LIFE OF THE MEDICAL DEVICE

**In order to limit negative impacts on the environment and to save natural resources, all waste must be treated according to its nature (recycling, recovery, incineration, landfill or other treatment for hazardous waste). In order to channel each waste to the appropriate treatment route, it is essential that it is collected and sorted correctly.**

The device must be dismantled, and the components must be handled according to the following requirements:

- Metal and plastic waste

Metal and plastic waste must be recycled according to the collection channels in place.

# ENZ

It must be recycled according to the collection channels in place.

They must be recycled and collected according to the channels in place.

For air dryers: Air dryers contain Hopcalite (which is a mixture of manganese dioxide and copper oxide; its current trade name is Carulite), aluminium oxide (or Alumina) and activated carbon. These products must be disposed of through a specific waste treatment process.

Therefore, the end-of-life management of electrical and electronic components is the total responsibility of the customer or the end-user.

The air unit with PROCOM2 or CYCLIC 2020 automate needs special precautions regarding EMC. The air unit with PROCOM2 or CYCLIC 2020 automate has to be installed and put into service according to the EMC requirements joined in the following tables.

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# ELECTROMAGNETIC COMPATIBILITY

The units with PROCOM2, PROCOM3, C2020 or MILLENIUM automate need special precautions regarding EMC.

## 1. Maintaining essential performance

Each source of vacuum, air or medical oxygen must maintain the pressure of the associated system within the required range. In the event of significant electromagnetic disturbances, the PLC controlling the control unit may malfunction. The worst cases are :

- continuous operation of the installation producing excess pressure or vacuum.
  - o In case of overpressure mechanical valves are present to protect the installation.
  - o Excess vacuum is not dangerous for patients.
- shutdown of the installation creating a lack of pressure or vacuum.
  - o For medical air, one or two backup mechanical pressure switches will operate the compressors.
  - o For medical vacuum, one or two backup mechanical vacuum switches will operate the vacuum pumps.
  - o For medical oxygen, a second source is present with similar equipment, then Mil's recommends the use of a 3rd source of gas supply operating by pressure difference, independently of any electronic system to ensure patient safety.

The maintenance indicated in the commissioning and maintenance manuals guarantees the integrity of the units and therefore their essential performances over time.

## 2. Instructions for use

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the unit, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

MIL'S ELECTROMEDICAL DEVICES are not compatible with HF ELECTROSURGICAL EQUIPMENT. As such, they must not be used in close proximity to each other.

## 3. Recommendations

The units with PROCOM2, PROCOM3, C2020 or MILLENIUM automate have to be installed and put into service according to the EMC requirements joined in the following tables.

Guidance and manufacturer's declaration – electromagnetic emissions – for all ME equipment and ME systems		
The unit with PROCOM2, PROCOM3, C2020 or MILLENIUM automate is intended for use in the electromagnetic environment specified below. The customer or the user of the unit with PROCOM2, PROCOM3, C2020 or MILLENIUM automate should assure that it is used in such an electromagnetic environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
Conducted emission CISPR 11 : 2009 / AMD1 : 2010	Group 1 / Class A COMPLIANT (Note A)	The units with PROCOM2, PROCOM3, C2020 or MILLENIUM automate uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Radiated emission CISPR 11 : 2009 / AMD1 : 2010	Group 1 / Class A COMPLIANT	
Harmonic emissions IEC 61000-3-2 : 2018	N/A	
Voltage fluctuations / flicker emissions IEC 61000-3-3 : 2013 / AMD1 : 2017	N/A	The units with PROCOM2, PROCOM3, C2020 or MILLENIUM automate is suitable for use in all establishments (e. g. hospitals), where the mains connection is isolated from the public low-voltage power supply network by transformers or substations.
Devices with motors or switching circuits CISPR 14-1	N/A	

# ELECTROMAGNETIC COMPATIBILITY

EN

Guidance and manufacturer's declaration – electromagnetic immunity - for all ME equipment and ME systems			
The unit with PROCOM2, PROCOM3, C2020 or MILLENIUM automate is intended for use in the electromagnetic environment specified below. The customer or the user of the unit with PROCOM2, PROCOM3, C2020 or MILLENIUM automate should assure that it is used in such an electromagnetic environment.			
Immunity test	IEC 60601 test level	Finding	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2 : 2008	±8 kV contact ±15 kV air	COMPLIANT (For MAVO, Note B and Note C)	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4 : 2012	±2 kV for power supply lines ±1 kV for input/output lines	COMPLIANT	Mains power quality should be that of a typical commercial or hospital environment.
Shock waves IEC 61000-4-5 : 2014 AMD1 : 2017	±1 kV mode differential ±2 kV mode common 0, 90, 180 & 270°	COMPLIANT (Note E)	
Transient electrical disturbances – 12 VDC power supply ISO 7637-2	600 V	N/A	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 : 2004 AMD1 : 2017	UT=0%, 0,5 cycle (0, 45, 90, 135, 180, 225, 270 & 315°) UT=0%; 1 cycle – UT=70% 25/30 cycles (0°) UT=0%; 250/300 cycles	COMPLIANT (Note H)	
Conducted RF IEC 61000-4-6 : 2013	3V (0.15-80MHz) 6V (discrete frequencies)	COMPLIANT (Note F)	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 : 2009	30 A/m 50 and 60 Hz	N/A (Note G)	No components sensible in product to magnetic field.
Radiated RF IEC 61000-4-3 : 2006 / AMD1 : 2007 / AMD2 : 2010	3 V/m 80 MHz to 2,7GHz  9V/m to 28 V/m 15 fréquences discrètes (tableau 9/ CEI 60601-1-2)	COMPLIANT (Note D)	<p>Portable and mobile RF communications equipment should be used no closer to any part of the units with PROCOM2, PROCOM3, C2020 or MILLENIUM automate, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p><math>d = 1.16 \sqrt{P}</math>  <math>d = 1.16 \sqrt{P}</math> from 80 MHz to 800 MHz  <math>d = 2.33 \sqrt{P}</math> from 800 MHz to 2,7 GHz  where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (1) should be less than the compliance level in each frequency range. (2)  Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
(1) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the units with PROCOM2, PROCOM3, C2020 or MILLENIUM automate is used exceeds the applicable RF compliance level above, the units with PROCOM2, PROCOM3, C2020 or MILLENIUM automate should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the units with PROCOM2, PROCOM3, C2020 or MILLENIUM automate.			
(2) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

# ELECTROMAGNETIC COMPATIBILITY

The units with PROCOM2, PROCOM3, C2020 or MILLENIUM automate can be affected by communication portables or mobiles RF's equipment.

Recommended separation distances between portable and mobile RF communications equipment and the unit with PROCOM2, PROCOM3, C2020 or MILLENIUM automate.			
The unit with PROCOM2, PROCOM3, C2020 or MILLENIUM automate is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the unit with PROCOM2, PROCOM3, C2020 or MILLENIUM automate can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit with PROCOM2, PROCOM3, C2020 or MILLENIUM automate as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz d = 1.17 √P	80 MHz to 800 MHz d = 1.17 √P	800 MHz to 2,7 GHz d = 2.33 √P
0,01	0.12 m	0.12 m	0.23 m
0,1	0.37 m	0.37 m	0.74 m
1	1.17 m	1.17 m	2.33 m
10	3.7 m	3.7 m	7.37 m
100	11.7 m	11.7 m	23.3 m
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Note A (see test report): Compliant with modifications either a Schaffner FN335-20-03 mains filter on each power supply (PROCOM and MILLENIUM) and a manufacturer filter on the inverter. Removing the original filter from the PROCOM power supply access.

Note B: On the Mavo report. Reset of the Mavo report then return to normal operation after a few seconds without operator intervention. Do not affect the essential criteria.

Note C: On the Mavo report, the electrostatic discharge (DES) on contact causes a reset of the Mavo report then a return to normal operation after a few seconds without operator intervention. Do not affect the essential criteria.

The ±15kV DES test on the Mavo report disabled it. However, the absence of the Mavo report does not affect the essential functions of the unit because this equipment is optional in the medical device.

On the flowmeter, the ESD test in contact causes the flowmeter to be out of order. However, this does not affect the essential functions of the unit because this equipment is optional in the medical device.

Note D: Discrete frequencies: 385MHz / 27V/m ; 450MHz / 28V/m ; 710MHz / 9V/m ; 745MHz / 9V/m ; 780MHz / 9V/m ; 810MHz / 28V/m ; 870MHz / 28V/m ; 930MHz / 28V/m ; 1720MHz / 28V/m ; 1845MHz / 28V/m ; 1970MHz / 28V/m ; 2450MHz / 28V/m ; 5240MHz / 9V/m ; 5500MHz / 9V/m ; 5785 MHz / 9V/m.

Note E: Last edition produced instead of the 2005 edition.

Note F: Discrete frequencies: 6.765MHz to 6.795MHz; 13.553MHz to 13.567MHz; 26.957MHz to 27.283MHz and 40.66MHz to 40.70MHz.

Note G: The equipment under test does not contain components sensitive to magnetic fields.

Note H: Deviation from the method. Considering the test means as well as the nature of the power supplies of the equipment under test, only the control elements are tested separately in single phase, i.e. PROCOM power supply, drive control and additional cabinet power supply.

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# SAFETY REGULATIONS

## MEANING OF THE PICTOGRAMS used in this manual and on machines

Pictograms have been inserted into the manual to make it clearer and to indicate particular points that should be observed or taken into account, as follows:

**Reading symbol: Mandatory:** Read the operator instructions (ISO 7000-M002)

**Mandatory:** Wear ear protectors (ISO 7010-M003)

**Warning:** Low temperature (ISO 7010 W010).

**Warning:** Maintenance operations are carried out (read the operator instructions).

**Warning:** High pressure in circuits and parts.

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**Warning:** Hot or dangerous gas exhaust point

**Warning:** The unit is remotely controlled and may start running without any warning being given (ISO 7000-W018).

**Warning:** risk of electric shock (ISO 7010-W012).

**Biological hazard symbol:** This pictogram informs of a risk or a danger, it is necessary to take care to handle the biological products or waste being able to infect the living organisms by inhalation, ingestion or cutaneous penetration.

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**Prohibition:** to remove the protective screen and the safety devices. (ISO 7010-M027)

**Warning:** do not touch the hot surfaces (risk of burns). (ISO 7010-W017)

Earth of the unit

Oil filling point (ISO 7000 0248)

**Environment protection symbol:** This pictogram calls to mind the requirement, during maintenance operations, to sort the waste, store it in a safe place and dispose of it in an environmentally-responsible manner.

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Lift-truck

Lifting point (ISO 7000 1368)

The removal of any self-adhesive label bearing the above pictograms from the machine is strictly forbidden. Any label which is damaged or peeling off must be replaced.

## 1. SAFETY REGULATIONS

**Would ask you to read the whole of this notice before commissioning your machine or carrying out any maintenance work on it. Our guarantee will not cover any damage caused by an unapproved siting or installation**

### 1.1. Warnings

The machinery described in this notice must not be used if there is any risk that the health of human beings, animals, plants and the environment may thereby be imperilled or destroyed.

The users and people responsible for servicing the machine must be fully responsible for the following:

- checking that the equipment is correctly installed and that good practice is followed when commissioning so that it is only used for its intended purpose,
- checking that the use of the equipment complies with regulations and that its operating conditions match the contractual specifications,
- ensuring compliance with local, regional and national regulations,
- ensuring compliance with the safety, operating and maintenance regulations laid down in this manual,
- ensuring that this manual is available to the staff servicing the equipment and ensuring that the staff observe the instructions, warnings and regulations set down in the manual.

Should the equipment be sent to the manufacturer for any reason, notification must be made indicating whether it is free of all substances harmful to health or to the contrary, that it is contaminated. If the equipment is contaminated, the contaminant and associated dangers must be declared.

Any machine arriving without a non-contamination statement, or statement indicating the nature of the contamination, will be sent back as is to the user. A full statement must accompany the machine. This contamination statement is mandatory to comply with legislation regarding health and safety. If you open a machine in your workshops, you are become accountable for any contamination which may be occasioned.

## 1.2. General informations

- This equipment is not suitable for use in areas where there is a risk of ignition or explosion.
- When shipping the equipment, account should be taken of local highway regulations if you have to use the public highway.
- Your safety and that of the people around you rely on the proper usage of the equipment. Ensure that you are completely familiar with the position and functions of all its controls.
- Any person using the machine or carrying out maintenance operations on it must be fully acquainted with it. Under no circumstances should any operation or other action be carried out unless the individual has the prerequisite skills.
- If you have to carry out unusual operations or work, contact your nearest dealer.
- He will put you in contact with specialists for this type of equipment.
- It is essential that all concerned know the location of the nearest fire extinguisher, emergency exits and first aid box to the machinery.
- Any person using the machinery or carrying out maintenance operations on it must be physically and mentally fit for the purpose. Said person must not be under the influence of alcohol, drugs or medication which could impair his reactions.
- Take the trouble to read and understand all the pictograms which indicate danger or precautions to take on the equipment.
- Avoid accidents which could be engendered by wearing clothing that is too loose or offering inadequate protection. Wear clothing that runs no risk of getting caught on parts. If necessary, wear additional protective gear (helmet, safety shoes, goggles, protective mask, gloves, reflective clothing, etc.)
- Because of the high decibel output of the machinery, the wearing of individual ear protection material is mandatory.
- As the machine may reach running temperatures in the region of 90°C (194°F), take all possible precautions against the risk of burns.
- The load bearing capacity of the floor on which the machinery is installed must be sufficient to bear its weight.
- If the floor is vulnerable to vibrations, shock absorbers and flexible pipes will have to be fitted to the machine.

## 1.3. General information on isolation and energy dissipation

Before any repair or servicing work is carried out, the machine must be fully or partially isolated and services to it cut.

Isolation and energy dissipation consists of:

- filling in the isolation log (summary of operations carried out) and informing the person responsible for isolation and energy dissipation (or the immediate superior) before taking the equipment out of service,
- putting the equipment into a specific state (cutting off the power, pneumatic, hydraulic and other energy supplies),
- keeping the equipment in this state (by padlocking, chaining, etc.),

- using panels to disclose the out-of service status of the equipment which record the following information:
  - ⇒ the date of isolation and energy dissipation,
  - ⇒ the equipment taken out of service,
  - ⇒ the reasons for removal from service.

The isolation and energy dissipation of machinery should only be reversed with the agreement of the person responsible for taking it out of service (or the immediate superior) and should be recorded in the isolating log.

## 1.4. Access and cleanliness

- Access to the machine must be unobstructed. Under no circumstances should an operator perform movements likely to cause personal injury to himself or damage property.
- The floor should be free of any trace or puddle of liquid of any description. Failure to observe this rule may lead to serious personal injuries.
- Oily rags and other inflammable items should be stored in a safe place.
- All environmentally responsible measures should be taken. In particular, all waste fluids should be recovered and recycled. Never dispose of these fluids directly into the environment.
- Call on the services of an accredited firm in the country where the equipment is used to recycle all oils, greases, and used filters.
- Fluids should never be poured into random containers, containers which are incorrectly labelled or open containers.
- Take care with all fluids used for servicing (oil, greases, etc.). They should not be drunk, inhaled or stored or conveyed in foodstuff containers (water bottles, beer bottles, etc.). Always use the original packing.

## 1.5. Safety checks (before starting up)

- For your safety and for the optimum service life of your machine, periodically inspect the machine from top to bottom before starting up.
- In particular, keep a watch out for loose bolts, missing bolts, build-up of dirt, spillage of liquids, etc.

## 1.6. Running


- If the machine is not running properly, or has parts missing or is out of-adjustment, stop the machine and repair it.  
Failure to observe this rule may lead to personal injuries.
- If a machine has not been fully repaired, it should not be re-started without the agreement of the immediate superior, who in giving his approval takes on all responsibility for incidents or accidents resulting.

## 1.7. Servicing-maintenance

- Always bear in mind that if servicing is not carried out carefully, observing all current recommendations, serious damage may occur which will be outside the terms of our guarantee.
- When servicing requires the work of several operators, there must be clear instructions for the tasks to be carried out. This is so that the manoeuvres and work between the members of the team are perfectly co-ordinated. Safety of all team members should be the prime concern of each individual team member.
- Never start repair work if you are unable to carry it out or complete it.
- During maintenance and servicing operations no unauthorised people should ever be permitted to work on the machinery.
- Always use the right tools for the tasks.

- Use cleaning products with care; they may be toxic and inflammable.
- The welding of parts onto the machine may damage mechanical and electrical elements. This sort of operation should not be carried out without consulting the manufacturer. Any disregard for this rule automatically nullifies the manufacturer's guarantee.
- For your safety and to avoid your own liability, only use MIL'S parts when replacing malfunctioning or broken down items.
- The immediate superior should be notified of all maintenance or testing operations on the equipment.
- During all maintenance operations the work area must be clearly delimited by white lines. At night, the area should be adequately lit.

- The lifting points are indicated by the  or  pictogram. During lifting operations, the machine must not be hoisted from any other point.
- If the machine has to be hoisted and moved for siting or maintenance operations, check that adequate lifting gear is used and familiarise yourselves with all the safety regulations surrounding the lifting of machinery.
- Under no circumstances should the machinery be needlessly suspended (during meal breaks, outside working hours, etc.)
- Following any work or repair on a machine function, carry out a check to see that this function has been rectified before using the machine.
- Extinguishers or firefighting equipment are the responsibility of the operator.

## 1.8. Planning maintenance operations

Many accidents can be avoided by organising maintenance operations. The supervisor must be clear about the tasks to be carried out and the possible risks involved. He should draw up a maintenance plan and explain it to all the people concerned so as to limit the risk of accidents.

The following pointers should be taken into account (non-exhaustive list):

- Which team members are needed for the work and what are their responsibilities?
- What is the weight of the load that may have to be lifted?
- What tooling is required for the job and how is it to be conveyed to the work area?
- What equipment is required to lift the equipment in absolute safety?

## 1.9. Electricity

- The supply voltage must not vary by more than 10%.
- An appropriate thermal relay is required for approved connection. The safety trip rating of the motor must equal the motor rating given on the manufacturer's plate.
- After connecting the motor or after each reconnection use a phase indicator to check the direction of rotation of the motor.
- The machine **must** be powered down whenever an electrical component is replaced.
- Before carrying out electrical connections on the machine, check that the voltage and supply frequency prevailing in the country of usage are identical to the operating voltage of the machine.
- The electricians commissioned to set up the machine are entirely responsible for the sizing of the electricity supply cable as well as all external cabling.
- The electrical connection of the machine must be carried out by a qualified electrician in compliance with regulations prevailing in the country where the machine is in use.

- Before carrying out any breakdown repairs or servicing, the machine must be partially or totally isolated and taken out of service.
- Any person carrying out servicing or breakdown repairs on a powered machine, must be qualified and suitably skilled to carry out such operations within that context.
- The machine must be partially taken out of service to prevent any other person starting it up.
- In the case of an electrical power cut or the tripping of a safety device, automatic restart of the machine must be disabled under all circumstances.
- Do not carry out any electrical connection or disconnection without knowledge of the circuitry. A faulty connection may cause personal injury or damage.
- Never wear any rings, bracelets or any other metallic object likely to earth a powered circuit, thereby causing burns.
- Following any work or repair on a machine function, check that this function has been rectified before using the machine.
- If maintenance operations are carried out on the machine while powered up, sparks may be produced on the electrical circuit. Do not work in confined areas in the presence of inflammable material, dust or vapour which could cause explosions or fires. Failure to observe this rule could result in personal or fatal injury.

## 1.10. Fluids

- A vacuum pump must be connected to suction piping without any mechanical stress. You are recommended to install flexible piping between the pump and the installation.
- A compressor, wherever possible, should be connected up for use by a pressure resistant flexible hose.
- The suction pressure of vacuum pumps must never exceed atmospheric pressure.
- The cross-section of vacuum pump exhaust lines must be at least as great as that of the pump connection hole. Too small a diameter on the line may create unacceptable pressure in the pump.
- No part of the casing should be exposed to vacuum or pressure.
- Do not start a vacuum pump running if its exhaust lines are blocked or constricted. Before starting the pump, check that the exhaust line is open and that it is free of any deposit that may have built up inside it.
- The choice of pipe materials should be appropriate for the process fluid.

Any serious incident occurring in connection with the device must be notified to the manufacturer and to the competent authority of the Member State in which the user is established.

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